

Dashboard for Monitoring and Signaling Elder's Heart Disease



Outline



01 Review Requirement

Review Requirement

Stakeholders

 Hospital, Organizations, and Families

Users

 Doctors , Organizations , Elders, and Families

Your goal

 To identify elder's heart disease before it happen

The input to your system

 Behavior data, sensor data and medical inspection result

The output of your system

Dashboard and Signal

Technologies do you need to achieve your goal

Big data and Machine learning

The goal of your team

What is the problem in Thailand you tried to solve?

 We choose Acute Coronary Syndrome-ACS in Elder.

How did you attempt to solve that problem?


 We create a tool to identify the problem before it happen.


What impact on the society did you expect your system to have?

 We expect that this tool can decrease deaths in Acute Coronary Syndrome-ACS more than 10%.

02 The System

Use case





Visualization

Query design

Heart disease risk

(Gauge)

BMI (Numeric) Body fat
percentage
(Donut)

Heart rate & Time (Line graph)

Differential
Pressure & Time
(Line graph)

& Time
(Line graph)

Dashboard design




Dashboard Link

https://app.powerbi.com/view?r=eyJrljoiMTlmZWYzZjAtMDE1Yi00NTM1LWlyYWlt MWRIMzliNGI0ODU4liwidCl6ljZmNDQzMmRjLTlwZDltNDQxZC1iMWRiLWFjMzM4 MGJhNjMzZClslmMiOjEwfQ%3D%3D&embedImagePlaceholder=true

03 Testing


Operation test result

Confusion Matrix



10-Fold Cross Validation

	Accuracy	AUC	Recall	Prec.	F1	Kappa	MCC
0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
1	1.0	1.0	1.0	1.0	1.0	1.0	1.0
2	1.0	1.0	1.0	1.0	1.0	1.0	1.0
3	1.0	1.0	1.0	1.0	1.0	1.0	1.0
4	1.0	1.0	1.0	1.0	1.0	1.0	1.0
5	1.0	1.0	1.0	1.0	1.0	1.0	1.0
6	1.0	1.0	1.0	1.0	1.0	1.0	1.0
7	1.0	1.0	1.0	1.0	1.0	1.0	1.0
8	1.0	1.0	1.0	1.0	1.0	1.0	1.0
9	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Mean	1.0	1.0	1.0	1.0	1.0	1.0	1.0
SD	0.0	0.0	0.0	0.0	0.0	0.0	0.0



Feature Importance

The way of the test

White Box Testing

Project Name: Dashboard for monitoring and signaling elder's heart disease Module Name: Pipeline and Dashboard Created by: Natchapol and Natchariya Date of creation: 1/3/2021 Date of review: 12/05/2021 to 20/05/2021



						<u>DR.FUZE</u>										
TEST CASE ID	TEST SCENARIO	TEST CASE	PRE-CONDITION	TEST STEPS	TEST DATA	EXPECTED RESULT	POST CONDITION	ACTUAL RESULT	STATUS (PASS/FAIL)							
			1. Access permission	1. Set trigger time		-: "										
1	Trigger pipeline	Trigger pipeline every 1 minute	to data factory	2. Click trigger		Pipeline run successfully after 1 minutes	Data in blob storage updated after 1 minutes	Pipeline run successfully after 1 minutes	Pass							
		Timilace		3. Wait for result shown		Timides	1 minutes									
		Trigger pipeline eveny	1. Access permission	1. Set trigger time		Dipoline rup guagagafully after	Data in blob storage undated after									
2	Trigger pipeline	Trigger pipeline every 5 minute	to data factory	2. Click trigger		Pipeline run successfully after 5 minutes	Data in blob storage updated after 5 minutes	Pipeline run successfully after 5 minutes	Pass							
		5 minute		3. Wait for result shown		3 minutes	5 Illinates									
		Trigger pipeline every	1. Access permission	1. Set trigger time		Pipeline run successfully after	Data in blob storage updated after									
3	Trigger pipeline	15 minute	to data factory	2. Click trigger		15 minutes	15 minutes	Pipeline run successfully after 15 minutes	Pass							
		15 mindec		3. Wait for result shown		15 mindees	15 mindees									
			Specify person ID Dackboard view of		Dashboard view change to											
4	Dashboard	Change person ID				that person ID		Dashboard view change to that person ID	Pass							
		-11.		Specify range of												
5	Dashboard	Filter range of		Datetime		Dashboard view change to		Dashboard view change to that range of datetime	Pass							
l		Datetime				that range of datetime		dateume								
			1. Choose person ID	1. Specify person ID			<u> </u>									
			2.Heart rate and	2. Wait for Heart rate and		1										
			Differential pressure	Differential pressure play												
_		Adjust value at	play finished	finished	Dashboard view Ad		Graph on Heart rate and Differential	Dashboard view Adjust value at Average	Page							
6	Dasilboard	Average of blood pressure		3. Click on Average of		at Average of blood pressure	pressure played finish.	of blood pressure	Pass							
		blood pressure graph														
			4. Adjust pressure value													
				5. Wait for result shown												





To ensure that there is no problem about connection.




What did you want to confirm through the test?

To ensure that all function work and connected.


04 Conclusion

01 Have you achieved your goal?





Achievement degree of your system




How did you Improve your system or development process when you faced problems?



What problem occurred?



Who and how solved that problem?



103 What are remaining issues for your system?

What is the remaining problem for your system?



How can you solve that problem in the future?



O5 Schedule and Role assignment

Finished

						Wor	king	Pro	ces	5										
Project Lead : Mr.Natchapo	oject Lead : Mr.Natchapol Patamawisut (Bank)													1/3/2	2021					Done
Secretary : Ms.Natchariya \					Finis	sh Pr	oject	: 31/	5/2021					In Process						
																	Late			
Scurm Master : Kao & Pro	rm Master : Kao & Project Manager : Bank																Not assign			
Work description	Assign to	Start Date	Mar-21					Ap	r-21	21		Ma	y-21		Percent	Finish Da	Period			Note
Work description	Assigit to	Glair Date	1	2	3	4	1	2	3	4	1	2	3	4	(%)	TIIISII De	116	Days		Note
Sprint1 : Conceptual design	n (1/3/2021 - 31	/3/2021)																		
Get requirement	All member	1/3/2021													100%	3/3/202	1	2		
Midterm Exam	All member	8/3/2021													NA	12/3/202	21	5	М	idterm Exam
Requirement Document	Kao,Bank	15/3/2021													100%	20/3/202	21	5		
Meeting	All member	27/3/2021													100%	28/3/202	21	1		
Presentation	All member	31/3/2021													100%	31/3/202	21	1	E	End Sprint 1

Finished

Sprint2 : System design (31/3/2021 - 26/4/2021)																				
Meeting	All member	31/3/2021										П	П	П	\neg		100%	31/3/2021	1	
																		I: 4/4/2021 (II:		
																		10/4/2021)(III :		Blob Storage ,
Create IaC	Eye	31/3/2021															100%	18/4/2021)(N:	~ 19	PostgreSQL and
																		21/4/2021)(V:		Data warehouse
												\dashv	\dashv	-	\dashv			25/4/2021		Manage data to
Programming	Bank	31/3/2021															100%	4/4/2021	5	Manage data to match the data
		041010004										\dashv	\dashv	\dashv	_		1000/	4442004		materiale data
Demo Visualization	Kao,New	31/3/2021										Ц					100%	4/4/2021	5	
Demo Docker	Fuze	31/3/2021															100%	5/4/2021	6	
Meeting	All member	4/4/2021															100%	4/4/2021	1	Update work
Extract and load Data	Fuze	5/4/2021															100%	21/4/2021	16	Upload & Download
Visualization	Kao,New	5/4/2021															100%	18/4/2021	13	Dashboard
Power BI POST request API	Као	12/4/021															100%	22/4/2021	10	
Machine Learning	Bank	5/4/2021															100%	18/4/2021	13	
Testing & Report	Kao,Bank	19/4/2021															100%	20/4/2021	1	
Presentation	All member	26/4/2021															100%	26/4/2021	1	End Sprint 2

Finished

Sprint3 : Implement & D	print3 : Implement & Deployment (26/4/2021 - 10/5/2021)																			
Meetings	All member	26/4/2021			Τ		Т	Т									100%	26/4/2021	1	
New Model	Bank	2/5/2021															100%	9/5/2021	7	
Airflow	Bank	26/4/2021			T		Τ										100%	2/5/2021	6	Automate workflow
Deploy Airflow	Fuze	2/5/2021															100%	9/5/2021	7	
Data Factories	Kao	2/5/2021															100%	11/5/2021	93	Transform Data
Data Label	Eye	2/5/2021			Ι		Ι										100%	7/5/2021	5	
Testing	Kao,Bank	4/5/2021			T												100%	11/5/2021	7	
End Sprint 3	All member	10/5/2021															100%	10/5/2021	1	End Sprint 3

Working Step

Sprint 4 : Testing (10/5/202	Sprint 4 : Testing (10/5/2021 - 7/6/2021)																			
Meetings	All member	12/5/2021															100%	12/5/2021	1	
Testing Model	Bank	12/5/2021															100%	20/5/2021	8	
Testing Pipeline	Eye	12/5/2021															100%	20/5/2021	8	
Testing Infrastructure	Fuze	12/5/2021															100%	20/5/2021	8	
Testing Dashboard	New	12/5/2021															100%	20/5/2021	8	
Testing & Report	Kao,Bank	19/5/2021															100%	21/5/2021	3	All testing & Report
Final Presentation	All member	5/6/2021															100%	5/6/2021	1	End Sprint 4

Members & Responsibility



Ms.Natchariya Wongamnuayporn (Kao) 61070507204



Mr.Natchapol Patamawisut (Bank) 61070507205



Ms.Thanaporn Cheenthada (Eye) 61070507209



Mr. Rungwigrai Payakkanuwat (New) 61070507219



Mr. Prakasit Chaiphet (Fuze) 61070507230

Responsibility

Bank: Leader, Programmer

Kao: Secretary, Tester

Eye: Infrastructure

New: Designer

Fuze: Deployment

Reference

- https://www.siphhospital.com/th/news/article/share/857/Acutecoronarysyndrome
- https://www.mayoclinic.org/diseases-conditions/acute-coronary-syndrome/symptoms-causes/syc-20352136#:~:text=Acute%20coronary%20syndrome%20is%20a,damaged%20or%20destroyed%20heart%20tissue.
- https://www.heart.org/en/health-topics/heart-attack/about-heart-attacks/acute-coronary-syndrome